

Product Datasheet

Human TK1 ELISA Kit (orb406721)

Catalog Number	orb406721
Category	Assays and Kits
Description	The Human TK1 ELISA Kit detects the target TK1 in human samples. It measures levels in sample types including serum, plasma and tissue homogenates. The detection range is 62.5 pg/mL-4000 pg/mL. The kit offers a sensitivity of 15.6 pg/mL.
Target	TK1
Reactivity	Human
UniProt ID	P04183
Application notes	Detection Wavelength: 450 nm
Assay Type	Sandwich
Assay Time	1-5h
Range	62.5 pg/mL-4000 pg/mL
Sensitivity	15.6 pg/mL
Sample Types	serum, plasma, tissue homogenates
Note	For research use only
Expiration Date	Please enquire.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)