

Product Datasheet

Human RF antibody (IgA) ELISA Kit (orb406748)

Catalog Number	orb406748
Category	Assays and Kits
Description	The Human RF antibody (IgA) ELISA Kit detects the target RF Ab (IgA) in human samples. Sample types including serum and plasma.
Target	RF Ab (IgA)
Reactivity	Human
Application notes	Detection Wavelength: 450 nm
Assay Type	Indirect
Assay Time	1-5h
Range	Request Information
Sensitivity	Request Information
Sample Types	serum, plasma
Note	For research use only
Expiration Date	Please enquire.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)