

Product Datasheet

Human ANXA2 ELISA Kit (orb406838)

Catalog Number	orb406838
Category	Assays and Kits
Description	The Human ANXA2 ELISA Kit detects the target ANXA2 in human samples. It measures levels in sample types including serum and urine and saliva and cell lysates. The detection range is 0.39 ng/mL-25 ng/mL. The kit offers a sensitivity of 0.098 ng/mL.
Target	ANXA2
Reactivity	Human
UniProt ID	P07355
Application notes	Detection Wavelength: 450 nm
Assay Type	Sandwich
Assay Time	1-5h
Range	0.39 ng/mL-25 ng/mL
Sensitivity	0.098 ng/mL
Sample Types	serum, urine, saliva, cell lysates
Note	For research use only
Expiration Date	Please enquire.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)