

Product Datasheet

Human ENK ELISA Kit (orb406854)

Catalog Number	orb406854
Category	Assays and Kits
Description	The Human ENK ELISA Kit detects the target ENK in human samples. It measures levels in sample types including serum and plasma. The detection range is 0.05 ng/mL-10 ng/mL. The kit offers a sensitivity of 0.025 ng/mL.
Target	ENK
Reactivity	Human
Application notes	Detection Wavelength: 450 nm
Assay Type	Competitive
Assay Time	1-5h
Range	0.05 ng/mL-10 ng/mL
Sensitivity	0.025 ng/mL
Sample Types	serum, plasma
Note	For research use only
Expiration Date	Please enquire.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)