

Product Datasheet

Human SERPINA6 ELISA Kit (orb406899)

Catalog Number	orb406899
Category	Assays and Kits
Description	The Human SERPINA6 ELISA Kit detects the target SERPINA6 in human samples. It measures levels in sample types including serum, plasma and tissue homogenates. The detection range is 20 pmol/mL-8000 pmol/mL. The kit offers a sensitivity of 20 pmol/mL.
Target	SERPINA6
Reactivity	Human
UniProt ID	P08185
Application notes	Detection Wavelength: 450 nm
Assay Type	Competitive
Assay Time	1-5h
Range	20 pmol/mL-8000 pmol/mL
Sensitivity	20 pmol/mL
Sample Types	serum, plasma, tissue homogenates
Note	For research use only
Expiration Date	Please enquire.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)