

Product Datasheet

Human PROK2 ELISA Kit (orb406916)

Catalog Number	orb406916
Category	Assays and Kits
Description	The Human PROK2 ELISA Kit detects the target PROK2 in human samples. It measures levels in sample types including serum, plasma, cell culture supernates and tissue homogenates. The detection range is 6.25 pg/mL-400 pg/mL. The kit offers a sensitivity of 1.56 pg/mL.
Target	PROK2
Reactivity	Human
UniProt ID	Q9HC23
Application notes	Detection Wavelength: 450 nm
Assay Type	Sandwich
Assay Time	1-5h
Range	6.25 pg/mL-400 pg/mL
Sensitivity	1.56 pg/mL
Sample Types	serum, plasma, cell culture supernates, tissue homogenates
Note	For research use only
Expiration Date	Please enquire.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)