

Product Datasheet

Human WARS ELISA Kit (orb407013)

Catalog Number	orb407013
Category	Assays and Kits
Description	The Human WARS ELISA Kit detects the target WARS in human samples. It measures levels in sample types including serum, plasma, cell lysates and tissue homogenates. The detection range is 0.312 ng/mL-20 ng/mL. The kit offers a sensitivity of 0.078 ng/mL.
Target	WARS
Reactivity	Human
UniProt ID	P23381
Application notes	Detection Wavelength: 450 nm
Assay Type	Sandwich
Assay Time	1-5h
Range	0.312 ng/mL-20 ng/mL
Sensitivity	0.078 ng/mL
Sample Types	serum, plasma, cell lysates, tissue homogenates
Note	For research use only
Expiration Date	Please enquire.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)