

Product Datasheet

Human CD106 ELISA Kit (orb407052)

Catalog Number	orb407052
Category	Assays and Kits
Description	The Human CD106 ELISA Kit detects the target L1CAM in human samples. It measures levels in sample types including serum, plasma and tissue homogenates. The detection range is 1.25 ng/mL-80 ng/mL. The kit offers a sensitivity of 0.31 ng/mL.
Target	L1CAM
Reactivity	Human
UniProt ID	P32004
Application notes	Detection Wavelength: 450 nm
Assay Type	Sandwich
Assay Time	1-5h
Range	1.25 ng/mL-80 ng/mL
Sensitivity	0.31 ng/mL
Sample Types	serum, plasma, tissue homogenates
Note	For research use only
Expiration Date	Please enquire.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)