

Product Datasheet

Human NOS1 ELISA Kit (orb407054)

Catalog Number	orb407054
Category	Assays and Kits
Description	The Human NOS1 ELISA Kit detects the target NOS1 in human samples. It measures levels in sample types including serum, plasma and cell culture supernates. The detection range is 0.156 IU/mL-10 IU/mL. The kit offers a sensitivity of 0.039 IU/mL.
Target	NOS1
Reactivity	Human
UniProt ID	P29475
Application notes	Detection Wavelength: 450 nm
Assay Type	Sandwich
Assay Time	1-5h
Range	0.156 IU/mL-10 IU/mL
Sensitivity	0.039 IU/mL
Sample Types	serum, plasma, cell culture supernates
Note	For research use only
Expiration Date	Please enquire.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)