

Product Datasheet

Human BRAF ELISA Kit (orb407101)

Catalog Number	orb407101
Category	Assays and Kits
Description	The Human BRAF ELISA Kit detects the target BRAF in human samples. It measures levels in sample types including serum and plasma and cell lysates. The detection range is 31.25 pg/mL-2000 pg/mL. The kit offers a sensitivity of 7.8 pg/mL.
Target	BRAF
Reactivity	Human
UniProt ID	P15056
Application notes	Detection Wavelength: 450 nm
Assay Type	Sandwich
Assay Time	1-5h
Range	31.25 pg/mL-2000 pg/mL
Sensitivity	7.8 pg/mL
Sample Types	serum, plasma, cell lysates
Note	For research use only
Expiration Date	Please enquire.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)