

Product Datasheet

Human PMP2 ELISA Kit (orb407121)

| | |
|--------------------------|--|
| Catalog Number | orb407121 |
| Category | Assays and Kits |
| Description | The Human PMP2 ELISA Kit detects the target PMP2 in human samples. Assay time: 1-5h. |
| Target | PMP2 |
| Reactivity | Human |
| UniProt ID | P02689 |
| Application notes | Detection Wavelength: 450 nm |
| Assay Time | 1-5h |
| Range | 0.312 ng/ml-20 ng/ml |
| Sensitivity | 0.078ng/ml |
| Sample Types | serum, plasma, tissue homogenates, cell lysates |
| Note | For research use only |
| Expiration Date | Please enquire. |

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)