

Product Datasheet

Human desmosine ELISA Kit (orb407132)

Catalog Number	orb407132
Category	Assays and Kits
Description	The Human desmosine ELISA Kit detects the target DES in human samples. It measures levels in sample types including serum and urine and tissue homogenates. The detection range is 0.156 ng/mL-10 ng/mL. The kit offers a sensitivity of 0.039 ng/mL.
Target	DES
Reactivity	Human
Application notes	Detection Wavelength: 450 nm
Assay Type	Sandwich
Assay Time	1-5h
Range	0.156 ng/mL-10 ng/mL
Sensitivity	0.039 ng/mL
Sample Types	serum, urine, tissue homogenates
Note	For research use only
Expiration Date	Please enquire.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)