

Product Datasheet

Human CSH1 ELISA Kit (orb407135)

Catalog Number	orb407135
Category	Assays and Kits
Description	The Human CSH1 ELISA Kit detects the target CSH1 in human samples. It measures levels in sample types including serum and plasma. The detection range is 80 ng/mL-4000 ng/mL. The kit offers a sensitivity of 100 ng/mL.
Target	CSH1
Reactivity	Human
UniProt ID	PODML2
Application notes	Detection Wavelength: 450 nm
Assay Type	Sandwich
Assay Time	1-5h
Range	80 ng/mL-4000 ng/mL
Sensitivity	100 ng/mL
Sample Types	serum, plasma
Note	For research use only
Expiration Date	Please enquire.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)