

Product Datasheet

Human ALPP ELISA Kit (orb407137)

Catalog Number	orb407137
Category	Assays and Kits
Description	The Human ALPP ELISA Kit detects the target ALPP in human samples. It measures levels in sample types including serum, plasma and tissue homogenates. The detection range is 84 pg/mL-2000 pg/mL. The kit offers a sensitivity of 34 pg/mL.
Target	ALPP
Reactivity	Human
UniProt ID	P05187
Application notes	Detection Wavelength: 450 nm
Assay Type	Sandwich
Assay Time	1-5h
Range	84 pg/mL-2000 pg/mL
Sensitivity	34 pg/mL
Sample Types	serum, plasma, tissue homogenates
Note	For research use only
Expiration Date	Please enquire.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)