

Product Datasheet

Human PTX3 ELISA Kit (orb407226)

Catalog Number	orb407226
Category	Assays and Kits
Description	The Human PTX3 ELISA Kit detects the target PTX3 in human samples. It measures levels in sample types including serum, plasma and tissue homogenates. The detection range is 78 pg/ml - 5000 pg/ml. The kit offers a sensitivity of 19.5 pg/ml.
Target	PTX3
Reactivity	Human
UniProt ID	P26022
Application notes	Detection Wavelength: 450 nm
Assay Type	Sandwich
Assay Time	1-5h
Range	78 pg/ml - 5000 pg/ml
Sensitivity	19.5 pg/ml
Sample Types	serum, plasma, tissue homogenates
Note	For research use only
Expiration Date	Please enquire.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)