

## Product Datasheet

### Human SEPP1 ELISA Kit (orb407228)

<b>Catalog Number</b>	orb407228
<b>Category</b>	Assays and Kits
<b>Description</b>	The Human SEPP1 ELISA Kit detects the target SEPP1 in human samples. It measures levels in sample types including serum, plasma, cell culture supernates and tissue homogenates and urine. The detection range is 0.156 ng/mL-10 ng/mL. The kit offers a sensitivity of 0.039 ng/mL.
<b>Target</b>	SEPP1
<b>Reactivity</b>	Human
<b>UniProt ID</b>	<b>P49908</b>
<b>Application notes</b>	Detection Wavelength: 450 nm
<b>Assay Type</b>	Sandwich
<b>Assay Time</b>	1-5h
<b>Range</b>	0.156 ng/mL-10 ng/mL
<b>Sensitivity</b>	0.039 ng/mL
<b>Sample Types</b>	serum, plasma, cell culture supernates, tissue homogenates, urine
<b>Note</b>	For research use only
<b>Expiration Date</b>	Please enquire.

**Biorbyt Ltd.**

7 Signet Court, Swann Road  
Cambridge  
CB5 8LA  
United Kingdom

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

**Biorbyt LLC**

68 TW Alexander Drive  
Research Triangle Park  
Durham  
NC 27713  
United States

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)