

Product Datasheet

Human CTLA4 ELISA Kit (orb407236)

Catalog Number	orb407236
Category	Assays and Kits
Description	The Human CTLA4 ELISA Kit detects the target CTLA4 in human samples. It measures levels in sample types including serum, plasma and tissue homogenates. The detection range is 125 pg/mL-8000 pg/mL. The kit offers a sensitivity of 31.25 pg/mL.
Target	CTLA4
Reactivity	Human
UniProt ID	P16410
Application notes	Detection Wavelength: 450 nm
Assay Type	Sandwich
Assay Time	1-5h
Range	125 pg/mL-8000 pg/mL
Sensitivity	31.25 pg/mL
Sample Types	serum, plasma, tissue homogenates
Note	For research use only
Expiration Date	Please enquire.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)