

## Product Datasheet

### Human KRT20 ELISA Kit (orb407244)

|                          |   |
|--------------------------|---|
| <b>Catalog Number</b>    | orb407244   |
| <b>Category</b>          | Assays and Kits   |
| <b>Description</b>       | The Human KRT20 ELISA Kit detects the target KRT20 in human samples. It measures levels in sample types including serum, plasma and tissue homogenates. The detection range is 0.625 ng/mL-40 ng/mL. The kit offers a sensitivity of 0.156 ng/mL. |
| <b>Target</b>            | KRT20   |
| <b>Reactivity</b>        | Human   |
| <b>UniProt ID</b>        | <b>P35900</b>   |
| <b>Application notes</b> | Detection Wavelength: 450 nm  |
| <b>Assay Type</b>        | Sandwich  |
| <b>Assay Time</b>        | 1-5h  |
| <b>Range</b>             | 0.625 ng/mL-40 ng/mL  |
| <b>Sensitivity</b>       | 0.156 ng/mL   |
| <b>Sample Types</b>      | serum, plasma, tissue homogenates   |
| <b>Note</b>              | For research use only   |
| <b>Expiration Date</b>   | Please enquire.   |

**Biorbyt Ltd.**

7 Signet Court, Swann Road  
Cambridge  
CB5 8LA  
United Kingdom

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

**Biorbyt LLC**

68 TW Alexander Drive  
Research Triangle Park  
Durham  
NC 27713-2847  
United States

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)