

Product Datasheet

Human CYP19A1 ELISA Kit (orb407252)

Catalog Number	orb407252
Category	Assays and Kits
Description	The Human CYP19A1 ELISA Kit detects the target CYP19A1 in human samples. This is a quantitative sandwich ELISA kit.
Target	CYP19A1
Reactivity	Human
UniProt ID	P11511
Application notes	Detection Wavelength: 450 nm
Assay Type	Sandwich
Assay Time	1-5h
Range	Request Information
Sensitivity	Request Information
Sample Types	Request Information
Note	For research use only

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)