

Product Datasheet

Human CYP3A4 ELISA Kit (orb407255)

Catalog Number	orb407255
Category	Assays and Kits
Description	The Human CYP3A4 ELISA Kit detects the target CYP3A4 in human samples. It measures levels in sample types including serum, plasma and tissue homogenates and cell lysates. The detection range is 23.5 pg/mL-1500 pg/mL. The kit offers a sensitivity of 5.8 pg/mL.
Target	CYP3A4
Reactivity	Human
UniProt ID	P08684
Application notes	Detection Wavelength: 450 nm
Assay Type	Sandwich
Assay Time	1-5h
Range	23.5 pg/mL-1500 pg/mL
Sensitivity	5.8 pg/mL
Sample Types	serum, plasma, tissue homogenates, cell lysates
Note	For research use only
Expiration Date	Please enquire.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)