

Product Datasheet

Human NPPA ELISA Kit (orb407295)

Catalog Number	orb407295
Category	Assays and Kits
Description	The Human NPPA ELISA Kit detects the target NPPA in human samples. It measures levels in sample types including serum and cell culture supernates. The detection range is 31.25 pg/mL-2000 pg/mL. The kit offers a sensitivity of 7.8 pg/mL.
Target	NPPA
Reactivity	Human
UniProt ID	P01160
Application notes	Detection Wavelength: 450 nm
Assay Type	Sandwich
Assay Time	1-5h
Range	31.25 pg/mL-2000 pg/mL
Sensitivity	7.8 pg/mL
Sample Types	serum, cell culture supernates
Note	For research use only

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)