

Product Datasheet

Human CXCL14 ELISA Kit (orb407310)

Catalog Number	orb407310
Category	Assays and Kits
Description	The Human CXCL14 ELISA Kit detects the target CXCL14 in human samples. It measures levels in sample types including serum, plasma and tissue homogenates. The detection range is 0.39 ng/ml-25 ng/ml. The kit offers a sensitivity of 0.09 ng/ml.
Target	CXCL14
Reactivity	Human
UniProt ID	O95715
Application notes	Detection Wavelength: 450 nm
Assay Type	Sandwich
Assay Time	1-5h
Range	0.39 ng/ml-25 ng/ml.
Sensitivity	0.09 ng/ml.
Sample Types	serum, plasma, tissue homogenates
Note	For research use only
Expiration Date	Please enquire.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)