

Product Datasheet

Human angiotensin II receptor 1 ELISA Kit (orb407324)

Catalog Number	orb407324
Category	Assays and Kits
Description	The Human angiotensin II receptor 1 ELISA Kit detects the target AT1R1 Ab in human samples. It measures levels in sample types including serum, plasma and tissue homogenates and cell lysates. The detection range is 0.312 ng/mL-20 ng/mL. The kit offers a sensitivity of 0.078 ng/mL.
Target	AT1R1 Ab
Reactivity	Human
Application notes	Detection Wavelength: 450 nm
Assay Type	Sandwich
Assay Time	1-5h
Range	0.312 ng/mL-20 ng/mL
Sensitivity	0.078 ng/mL
Sample Types	serum, plasma, tissue homogenates, cell lysates
Note	For research use only
Expiration Date	Please enquire.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)