

Product Datasheet

Human SERPING1 ELISA Kit (orb407355)

Catalog Number	orb407355
Category	Assays and Kits
Description	The Human SERPING1 ELISA Kit detects the target SERPING1 in human samples. It measures levels in sample types including serum, plasma and tissue homogenates. The detection range is 18.75 ng/mL-1200 ng/mL. The kit offers a sensitivity of 4.7 ng/mL.
Target	SERPING1
Reactivity	Human
UniProt ID	P05155
Application notes	Detection Wavelength: 450 nm
Assay Type	Competitive
Assay Time	1-5h
Range	18.75 ng/mL-1200 ng/mL
Sensitivity	4.7 ng/mL
Sample Types	serum, plasma, tissue homogenates
Note	For research use only
Expiration Date	Please enquire.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)