

Product Datasheet

Human SAA4 ELISA Kit (orb407360)

Catalog Number	orb407360
Category	Assays and Kits
Description	The Human SAA4 ELISA Kit detects the target SAA4 in human samples. It measures levels in sample types including serum, plasma and tissue homogenates. The detection range is 15.6 ng/ml-1000 ng/ml. The kit offers a sensitivity of 3.9 ng/ml.
Target	SAA4
Reactivity	Human
UniProt ID	P35542
Application notes	• Sample Type: serum, plasma, tissue homogenates • Assay Time: 1-5 hours • Sample Volume: 50-100 µl • Detection Wavelength: 450 nm
Assay Type	Sandwich
Assay Time	1-5h
Range	15.6 ng/ml-1000 ng/ml.
Sensitivity	3.9 ng/ml.
Sample Types	serum, plasma, tissue homogenates
Note	For research use only
Expiration Date	Please enquire.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)