

Product Datasheet

Human HPSE ELISA Kit (orb407438)

Catalog Number	orb407438
Category	Assays and Kits
Description	The Human HPSE ELISA Kit detects the target HPSE in human samples. It measures levels in sample types including serum and plasma and ascitic fluid and cell culture supernates and tissue homogenates. The detection range is 0.156 ng/mL-10 ng/mL. The kit offers a sensitivity of 0.039 ng/mL.
Target	HPSE
Reactivity	Human
UniProt ID	Q9Y251
Application notes	Detection Wavelength: 450 nm
Assay Type	Sandwich
Assay Time	1-5h
Range	0.156 ng/mL-10 ng/mL
Sensitivity	0.039 ng/mL
Sample Types	serum, plasma, ascitic fluid, cell culture supernates, tissue homogenates
Note	For research use only
Expiration Date	Please enquire.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)