

Product Datasheet

Human FBN1 ELISA Kit (orb407476)

Catalog Number	orb407476
Category	Assays and Kits
Description	The Human FBN1 ELISA Kit detects the target FBN1 in human samples. It measures levels in sample types including serum, plasma and tissue homogenates. The detection range is 0.235 ng/ml-15 ng/ml. The kit offers a sensitivity of 0.058 ng/ml.
Target	FBN1
Reactivity	Human
UniProt ID	P35555
Application notes	Detection Wavelength: 450 nm
Assay Type	Sandwich
Assay Time	1-5h
Range	0.235 ng/ml-15 ng/ml.
Sensitivity	0.058 ng/ml.
Sample Types	serum, plasma, tissue homogenates.
Note	For research use only
Expiration Date	Please enquire.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)