

Product Datasheet

Human APOH ELISA Kit (orb407488)

Catalog Number	orb407488
Category	Assays and Kits
Description	The Human APOH ELISA Kit detects the target APOH in human samples. It measures levels in sample types including serum, plasma and tissue homogenates. The detection range is 15.6 ng/mL-1000 ng/mL. The kit offers a sensitivity of 3.9 ng/mL.
Target	APOH
Reactivity	Human
UniProt ID	P02749
Application notes	Detection Wavelength: 450 nm
Assay Type	Sandwich
Assay Time	1-5h
Range	15.6 ng/mL-1000 ng/mL
Sensitivity	3.9 ng/mL
Sample Types	serum, plasma, tissue homogenates
Note	For research use only
Expiration Date	Please enquire.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)