

Product Datasheet

Human MUC2 ELISA Kit (orb407506)

Catalog Number	orb407506
Category	Assays and Kits
Description	The Human MUC2 ELISA Kit detects the target MUC2 in human samples. It measures levels in sample types including serum, plasma and cell culture supernates. The detection range is 3.12 ng/mL-200 ng/mL. The kit offers a sensitivity of 0.78 ng/mL.
Target	MUC2
Reactivity	Human
UniProt ID	Q02817
Application notes	Detection Wavelength: 450 nm
Assay Type	Sandwich
Assay Time	1-5h
Range	3.12 ng/mL-200 ng/mL
Sensitivity	0.78 ng/mL
Sample Types	serum, plasma, cell culture supernates
Note	For research use only
Expiration Date	Please enquire.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)