

Product Datasheet

Human MUC5AC ELISA Kit (orb407510)

Catalog Number	orb407510
Category	Assays and Kits
Description	The Human MUC5AC ELISA Kit detects the target MUC5AC in human samples. It measures levels in sample types including serum, plasma and cell culture supernates and urine. The detection range is 15.6 ng/mL-1000 ng/mL. The kit offers a sensitivity of 3.9 ng/mL.
Target	MUC5AC
Reactivity	Human
UniProt ID	P98088
Application notes	Detection Wavelength: 450 nm
Assay Type	Sandwich
Assay Time	1-5h
Range	15.6 ng/mL-1000 ng/mL
Sensitivity	3.9 ng/mL
Sample Types	serum, plasma, cell culture supernates, urine
Note	For research use only
Expiration Date	Please enquire.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)