

Product Datasheet

Human NEP ELISA Kit (orb407539)

Catalog Number	orb407539
Category	Assays and Kits
Description	The Human NEP ELISA Kit detects the target NEP in human samples. It measures levels in sample types including serum, plasma and tissue homogenates. The detection range is 78 pg/mL-5000 pg/mL. The kit offers a sensitivity of 19.5 pg/mL.
Target	NEP
Reactivity	Human
UniProt ID	P08473
Application notes	Detection Wavelength: 450 nm
Assay Type	Sandwich
Assay Time	1-5h
Range	78 pg/mL-5000 pg/mL
Sensitivity	19.5 pg/mL
Sample Types	serum, plasma, tissue homogenates
Note	For research use only

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)