

Product Datasheet

Human ADAM17 ELISA Kit (orb407544)

Catalog Number	orb407544
Category	Assays and Kits
Description	The Human ADAM17 ELISA Kit detects the target ADAM17 in human samples. It measures levels in sample types including serum, plasma and tissue homogenates and cell lysates. The detection range is 62.5 pg/mL-4000 pg/mL. The kit offers a sensitivity of 15.6 pg/mL.
Target	ADAM17
Reactivity	Human
UniProt ID	P78536
Application notes	Detection Wavelength: 450 nm
Assay Type	Sandwich
Assay Time	1-5h
Range	62.5 pg/mL-4000 pg/mL
Sensitivity	15.6 pg/mL
Sample Types	serum, plasma, tissue homogenates, cell lysates
Note	For research use only
Expiration Date	Please enquire.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)