

Product Datasheet

Human ANGPTL8 ELISA Kit (orb407621)

Catalog Number	orb407621
Category	Assays and Kits
Description	The Human ANGPTL8 ELISA Kit detects the target ANGPTL8 in human samples. It measures levels in sample types including serum, plasma and tissue homogenates. The detection range is 6.25 pg/mL-400 pg/mL. The kit offers a sensitivity of 1.56 pg/mL.
Target	ANGPTL8
Reactivity	Human
UniProt ID	Q6UXH0
Application notes	Detection Wavelength: 450 nm
Assay Type	Sandwich
Assay Time	1-5h
Range	6.25 pg/mL-400 pg/mL
Sensitivity	1.56 pg/mL
Sample Types	serum, plasma, tissue homogenates
Note	For research use only
Expiration Date	Please enquire.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)