

Product Datasheet

Mouse CLCF1 ELISA Kit (orb407844)

Catalog Number	orb407844
Category	Assays and Kits
Description	The Mouse CLCF1 ELISA Kit detects the target CLCF1 in mouse samples. It measures levels in sample types including serum, plasma, cell culture supernates and tissue homogenates. The detection range is 23.5 pg/mL-1500 pg/mL. The kit offers a sensitivity of 5.8 pg/mL.
Target	CLCF1
Reactivity	Mouse
UniProt ID	Q9QZM3
Application notes	Detection Wavelength: 450 nm
Assay Type	Sandwich
Assay Time	1-5h
Range	23.5 pg/mL-1500 pg/mL
Sensitivity	5.8 pg/mL
Sample Types	serum, plasma, cell culture supernates, tissue homogenates
Note	For research use only
Expiration Date	Please enquire.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)