

Product Datasheet

Mouse EGFL7 ELISA Kit (orb408057)

Catalog Number	orb408057
Category	Assays and Kits
Description	The Mouse EGFL7 ELISA Kit detects the target EGFL7 in mouse samples. It measures levels in sample types including serum and plasma and tissue homogenates. The detection range is 7.8 ng/mL-500 ng/mL. The kit offers a sensitivity of 1.95 ng/mL.
Target	EGFL7
Reactivity	Mouse
UniProt ID	Q9QXT5
Application notes	Detection Wavelength: 450 nm
Assay Type	Sandwich
Assay Time	1-5h
Range	7.8 ng/mL-500 ng/mL
Sensitivity	1.95 ng/mL
Sample Types	serum, plasma,tissue homogenates
Note	For research use only
Expiration Date	Please enquire.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)