

Product Datasheet

Mouse RSPO4 ELISA Kit (orb408592)

Catalog Number	orb408592
Category	Assays and Kits
Description	The Mouse RSPO4 ELISA Kit detects the target RSPO4 in mouse samples. It measures levels in sample types including serum, plasma and tissue homogenates. The detection range is 34.4 pg/mL-2200 pg/mL. The kit offers a sensitivity of 8.6 pg/mL.
Target	RSPO4
Reactivity	Mouse
UniProt ID	Q8BJ73
Application notes	Detection Wavelength: 450 nm
Assay Type	Sandwich
Assay Time	1-5h
Range	34.4 pg/mL-2200 pg/mL
Sensitivity	8.6 pg/mL
Sample Types	serum, plasma, tissue homogenates
Note	For research use only
Expiration Date	Please enquire.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)