

## Product Datasheet

### Mouse SFRP4 ELISA Kit (orb408601)

<b>Catalog Number</b>	orb408601
<b>Category</b>	Assays and Kits
<b>Description</b>	The Mouse SFRP4 ELISA Kit detects the target SFRP4 in mouse samples. It measures levels in sample types including serum, plasma, cell culture supernates and tissue homogenates. The detection range is 0.47 ng/mL-30 ng/mL. The kit offers a sensitivity of 0.11 ng/mL.
<b>Target</b>	SFRP4
<b>Reactivity</b>	Mouse
<b>UniProt ID</b>	<b>Q9Z1N6</b>
<b>Application notes</b>	Detection Wavelength: 450 nm
<b>Assay Type</b>	Sandwich
<b>Assay Time</b>	1-5h
<b>Range</b>	0.47 ng/mL-30 ng/mL
<b>Sensitivity</b>	0.11 ng/mL
<b>Sample Types</b>	serum, plasma, cell culture supernates, tissue homogenates
<b>Note</b>	For research use only
<b>Expiration Date</b>	Please enquire.

**Biorbyt Ltd.**

7 Signet Court, Swann Road  
Cambridge  
CB5 8LA  
United Kingdom

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

**Biorbyt LLC**

68 TW Alexander Drive  
Research Triangle Park  
Durham  
NC 27713  
United States

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)