

Product Datasheet

Mouse Interferon gamma ELISA Kit (orb408817)

Catalog Number	orb408817
Category	Assays and Kits
Description	The Mouse Interferon gamma ELISA Kit detects the target CXCL11 in mouse samples. It measures levels in sample types including serum, plasma and cell culture supernates and tissue homogenates. The detection range is 31.25 pg/mL-2000 pg/mL. The kit offers a sensitivity of 7.81 pg/mL.
Target	CXCL11
Reactivity	Mouse
UniProt ID	Q9JHH5
Application notes	Detection Wavelength: 450 nm
Assay Type	Sandwich
Assay Time	1-5h
Range	31.25 pg/mL-2000 pg/mL
Sensitivity	7.81 pg/mL
Sample Types	serum, plasma, cell culture supernates,tissue homogenates
Note	For research use only
Expiration Date	Please enquire.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)