

## Product Datasheet

### Mouse NAG ELISA Kit (orb408851)

<b>Catalog Number</b>	orb408851
<b>Category</b>	Assays and Kits
<b>Description</b>	The Mouse NAG ELISA Kit detects the target NAG in mouse samples. It measures levels in sample types including serum, plasma and tissue homogenates. The detection range is 15.6 $\mu$ U/mL-1000 $\mu$ U/mL. The kit offers a sensitivity of 3.9 $\mu$ U/mL.
<b>Target</b>	NAG
<b>Reactivity</b>	Mouse
<b>UniProt ID</b>	<b>P20060</b>
<b>Application notes</b>	Detection Wavelength: 450 nm
<b>Assay Type</b>	Sandwich
<b>Assay Time</b>	1-5h
<b>Range</b>	15.6 $\mu$ U/mL-1000 $\mu$ U/mL
<b>Sensitivity</b>	3.9 $\mu$ U/mL
<b>Sample Types</b>	serum, plasma, tissue homogenates
<b>Note</b>	For research use only
<b>Expiration Date</b>	Please enquire.

**Biorbyt Ltd.**

7 Signet Court, Swann Road  
Cambridge  
CB5 8LA  
United Kingdom

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

**Biorbyt LLC**

68 TW Alexander Drive  
Research Triangle Park  
Durham  
NC 27713  
United States

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)