

Product Datasheet

Mouse alpha 2 Macroglobulin ELISA Kit (orb408883)

| | |
|--------------------------|---|
| Catalog Number | orb408883 |
| Category | Assays and Kits |
| Description | The Mouse alpha 2 Macroglobulin ELISA Kit detects the target A2M in mouse samples. It measures levels in sample types including serum and plasma. The detection range is 0.156 µg/mL-10 µg/mL. The kit offers a sensitivity of 0.039 µg/mL. |
| Target | A2M |
| Reactivity | Mouse |
| UniProt ID | Q61838 |
| Application notes | Detection Wavelength: 450 nm |
| Assay Type | Competitive |
| Assay Time | 1-5h |
| Range | 0.156 µg/mL-10 µg/mL |
| Sensitivity | 0.039 µg/mL |
| Sample Types | serum, plasma |
| Note | For research use only |
| Expiration Date | Please enquire. |

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)