

## Product Datasheet

### Mouse MSH Alpha ELISA Kit (orb408908)

<b>Catalog Number</b>	orb408908
<b>Category</b>	Assays and Kits
<b>Description</b>	The Mouse MSH Alpha ELISA Kit detects the target POMC in mouse samples. It measures levels in sample types including serum, plasma, cell culture supernates and tissue homogenates. The detection range is 46.9 pg/ml-3000 pg/ml. The kit offers a sensitivity of 24 pg/ml.
<b>Target</b>	POMC
<b>Reactivity</b>	Mouse
<b>UniProt ID</b>	<b>P01193</b>
<b>Application notes</b>	Detection Wavelength: 450 nm
<b>Assay Type</b>	Competitive
<b>Assay Time</b>	1-5h
<b>Range</b>	46.9 pg/ml-3000 pg/ml.
<b>Sensitivity</b>	24 pg/ml
<b>Sample Types</b>	serum, plasma, cell culture supernates, tissue homogenates.
<b>Note</b>	For research use only
<b>Expiration Date</b>	Please enquire.

**Biorbyt Ltd.**

7 Signet Court, Swann Road  
Cambridge  
CB5 8LA  
United Kingdom

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

**Biorbyt LLC**

68 TW Alexander Drive  
Research Triangle Park  
Durham  
NC 27713-2847  
United States

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)