

Product Datasheet

Mouse TRH ELISA Kit (orb409001)

Catalog Number	orb409001
Category	Assays and Kits
Description	The Mouse TRH ELISA Kit detects the target TRH in mouse samples. It measures levels in sample types including serum. The detection range is 2.5 pg/mL-160 pg/mL. The kit offers a sensitivity of 0.625 pg/mL.
Target	TRH
Reactivity	Mouse
UniProt ID	Q62361
Application notes	Detection Wavelength: 450 nm
Assay Type	Sandwich
Assay Time	1-5h
Range	2.5 pg/mL-160 pg/mL
Sensitivity	0.625 pg/mL
Sample Types	serum
Note	For research use only
Expiration Date	Please enquire.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)