

Product Datasheet

mouse ACTH ELISA Kit (orb409003)

Catalog Number	orb409003
Category	Assays and Kits
Description	The mouse ACTH ELISA Kit detects the target ACTH in mouse samples. It measures levels in sample types including serum and plasma. The detection range is 18.75 pg/mL-1200 pg/mL. The kit offers a sensitivity of 4.68 pg/mL.
Target	ACTH
Reactivity	Mouse
Application notes	Detection Wavelength: 450 nm
Assay Type	Sandwich
Assay Time	1-5h
Range	18.75 pg/mL-1200 pg/mL
Sensitivity	4.68 pg/mL
Sample Types	serum, plasma
Note	For research use only
Expiration Date	Please enquire.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)