

Product Datasheet

Mouse FASLG ELISA Kit (orb409027)

Catalog Number	orb409027
Category	Assays and Kits
Description	The Mouse FASLG ELISA Kit detects the target FASLG in mouse samples. It measures levels in sample types including serum, plasma and tissue homogenates. The detection range is 7.8 pg/mL-500 pg/mL. The kit offers a sensitivity of 1.95 pg/mL.
Target	FASLG
Reactivity	Mouse
UniProt ID	P41047
Application notes	Detection Wavelength: 450 nm
Assay Type	Sandwich
Assay Time	1-5h
Range	7.8 pg/mL-500 pg/mL
Sensitivity	1.95 pg/mL
Sample Types	serum, plasma, tissue homogenates
Note	For research use only
Expiration Date	Please enquire.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)