

Product Datasheet

Mouse Cox2 ELISA Kit (orb409114)

Catalog Number	orb409114
Category	Assays and Kits
Description	The Mouse Cox2 ELISA Kit detects the target PTGS2 in mouse samples. It measures levels in sample types including serum, plasma and tissue homogenates and cell lysates. The detection range is 31.25 pg/mL-2000 pg/mL. The kit offers a sensitivity of 7.8 pg/mL.
Target	PTGS2
Reactivity	Mouse
UniProt ID	Q05769
Application notes	Detection Wavelength: 450 nm
Assay Type	Sandwich
Assay Time	1-5h
Range	31.25 pg/mL-2000 pg/mL
Sensitivity	7.8 pg/mL
Sample Types	serum, plasma, tissue homogenates, cell lysates
Note	For research use only
Expiration Date	Please enquire.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)