

Product Datasheet

Mouse CNTF ELISA Kit (orb409167)

Catalog Number	orb409167
Category	Assays and Kits
Description	The Mouse CNTF ELISA Kit detects the target CNTF in mouse samples. It measures levels in sample types including serum, plasma and tissue homogenates and cell culture supernates. The detection range is 6.25 pg/mL-400 pg/mL. The kit offers a sensitivity of 1.56 pg/mL.
Target	CNTF
Reactivity	Mouse
UniProt ID	P51642
Application notes	Detection Wavelength: 450 nm
Assay Type	Sandwich
Assay Time	1-5h
Range	6.25 pg/mL-400 pg/mL
Sensitivity	1.56 pg/mL
Sample Types	serum, plasma, tissue homogenates, cell culture supernates
Note	For research use only

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)