

Product Datasheet

Mouse MT-2 ELISA Kit (orb409171)

Catalog Number	orb409171
Category	Assays and Kits
Description	The Mouse MT-2 ELISA Kit detects the target MT-2 in mouse samples. It measures levels in sample types including serum, plasma and tissue homogenates. The detection range is 15.6 pg/mL-1000 pg/mL. The kit offers a sensitivity of 3.9 pg/mL.
Target	MT-2
Reactivity	Mouse
UniProt ID	P02798
Application notes	Detection Wavelength: 450 nm
Assay Type	Sandwich
Assay Time	1-5h
Range	15.6 pg/mL-1000 pg/mL
Sensitivity	3.9 pg/mL
Sample Types	serum, plasma, tissue homogenates
Note	For research use only
Expiration Date	Please enquire.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)