

Product Datasheet

Mouse Ceruloplasmin ELISA Kit (orb409405)

Catalog Number	orb409405
Category	Assays and Kits
Description	The Mouse Ceruloplasmin ELISA Kit detects the target CP in mouse samples. It measures levels in sample types including serum, plasma, cell culture supernates and tissue homogenates. The detection range is 31.2 ng/mL-2000 ng/mL. The kit offers a sensitivity of 7.8 ng/mL.
Target	CP
Reactivity	Mouse
UniProt ID	Q61147
Application notes	Detection Wavelength: 450 nm
Assay Type	Competitive
Assay Time	1-5h
Range	31.2 ng/mL-2000 ng/mL
Sensitivity	7.8 ng/mL
Sample Types	serum, plasma, cell culture supernates, tissue homogenates
Note	For research use only
Expiration Date	Please enquire.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)