

Product Datasheet

Mouse SOCS-1 ELISA Kit (orb409430)

Catalog Number	orb409430
Category	Assays and Kits
Description	The Mouse SOCS-1 ELISA Kit detects the target SOCS1 in mouse samples. It measures levels in sample types including serum, plasma and tissue homogenates. The detection range is 62.5 pg/mL-4000 pg/mL. The kit offers a sensitivity of 15.6 pg/mL.
Target	SOCS1
Reactivity	Mouse
UniProt ID	O35716
Application notes	Detection Wavelength: 450 nm
Assay Type	Sandwich
Assay Time	1-5h
Range	62.5 pg/mL-4000 pg/mL
Sensitivity	15.6 pg/mL
Sample Types	serum, plasma, tissue homogenates
Note	For research use only

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)