

Product Datasheet

Mouse SAA1 ELISA Kit (orb409468)

Catalog Number	orb409468
Category	Assays and Kits
Description	The Mouse SAA1 ELISA Kit detects the target SAA1 in mouse samples. It measures levels in sample types including serum and plasma and tissue homogenates. The detection range is 3.12 ng/ml - 200 ng/ml. The kit offers a sensitivity of 0.78 ng/ml.
Target	SAA1
Reactivity	Mouse
UniProt ID	P05366
Application notes	• Sample Type: serum, plasma, cell culture supernates, tissue homogenates• Assay Time: 1-5 hours• Sample Volume: 50-100 µl• Detection Wavelength: 450 nm
Assay Type	Sandwich
Assay Time	1-5h
Range	3.12 ng/ml - 200 ng/ml
Sensitivity	0.78 ng/ml
Sample Types	serum, plasma,tissue homogenates
Note	For research use only
Expiration Date	Please enquire.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)